TMX[®]
BEBO
BECKENBODENTRIGGER

TMX[®] **BEBO**

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TMX[®]
BEBO
PELVIC FLOOR TRIGGER

TMX[®] **BEBO**

OUR MOST EFFECTIVE TRIGGER
AGAINST PELVIC FLOOR PAIN AND DISCOMFORT

**WE MUST RELEASE BEFORE WE
STRENGTHEN!**

VERY EASY APPLICATION THAT DOES NOT DIFFER FOR MEN
AND WOMEN.

Let's not kid ourselves, it's a global issue. It can be assumed that pelvic floor complaints affect up to 90% women and the remaining 10% men.

A woman's pelvic floor suffers major trauma during childbirth. The pelvic floor muscles are stretched to the maximum and lose their functionality as a result. The physiological tension decreases and the controlled tensing and relaxing of the pelvic floor is severely impaired. Functional disorders, slight tearing and everyday activities become uncontrolled strains. In men, on the other hand, the pelvic floor muscles are most affected by inflammatory processes of the prostate or bladder. Age-related loss of tension and functional limitations are unfortunately an annoying side effect for both sexes.

The pelvic floor is nothing more than a combination of a large number of muscles that are understood as the diaphragm of the small pelvis. The functionality of the pelvic floor is essential for our everyday life, e.g. for controlled urination, bowel movements and also for sexual needs.

**APPLICATION
WITH THE TMX[®]
3-Step-METHOD:**

- 1. Triggering:** Place the TMX[®] BEBO on a chair and sit on it. The targeted pressure on your muscle allows you to counteract deep-seated tension.
- 2. Mobilize during triggering:** The subsequent mobilization gives the muscle back its freedom of movement. To do this, move your upper body slightly forwards and backwards.
- 3. Awareness:** Awareness of the pelvic floor muscles rounds off the relaxation process.

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


LOW INTENSITY

(between pubic bone and anus)
Trigger ridges to the front

HIGH INTENSITY

(between pubic bone and anus)
Large trigger to the front

TRIGGER REAR PERINEUM

LOW INTENSITY

(between anus and coccyx)
Trigger ridges to the front


TRIGGER REAR PERINEUM

HIGH INTENSITY


(between anus and coccyx)
Large trigger to the front

AREAS OF APPLICATION

FEMALE PELVIC FLOOR

MALE PELVIC FLOOR

1 ANTERIOR PERINEUM 2 REAR PERINEUM

INDICATIONS

Pelvic floor complaints, pelvic floor pain, pelvic floor tension, pelvic floor insufficiency, urination disorders, prostatitis, stress incontinence, muscular spasms, vaginismus

CONTRAINDICATIONS

Inflammation in the pelvic floor, signs of inflammation (warmth, swelling, redness), after operations, open wounds, cysts, scars

FREE TMX[®] WORLD OF EXERCISE

With the free TMX EXERCISE WORLD you get a variety of exercises for your pelvic floor.

Simply scan the QR code or start directly on uebungswelt.de!

Direct help from our physio team:
PHYSIO@TMX-TRIGGER.DE

TMXTRIGGER